

## FEATURES

Extensive range of options available



### ALSO AVAILABLE

- Frame widths up to 36" (91 cm)
- Custom frame depths
- Custom seats

SPECIFICATIONS	Raz-SP600 (Self Propelled 600 lbs capacity)
WEIGHT CAPACITY	600 lbs (275 kg)
FRAME WIDTH (seat and back frame)	22" (56 cm) [Custom widths available]
OVERALL WIDTH (with standard armpad)	30" (76 cm) with Standard Arm Pad
SEAT WIDTH (commode seat)	23" (58 cm)
SEAT DEPTH (front of back post to front of seat)	18" (46 cm)
OVERALL LENGTH WHEN UPRIGHT With Z709 Footrests	40" – 42" (101.5 – 106.5 cm)
CLEARANCE HEIGHT UNDER SEAT	5" (125 mm) casters with pan brackets: 15.25" or 16.25" (39 cm or 41 cm) 5" (125 mm) casters without pan brackets: 16" or 17" (40.5 cm or 43 cm) 6" (125 mm) casters with pan brackets: 16.25", 17.25", 20.25" (41 cm, 44 cm, 51.5 cm) 6" (125 mm) casters without pan brackets: 17", 18", or 21" (43 cm, 46 cm, 53 cm)
ARMREST HEIGHT ABOVE SEAT	9.5" (24 cm)
SEAT-TO-FLOOR HEIGHT (5" casters)	20" and 21" (51 cm and 53 cm)
SEAT-TO-FLOOR HEIGHT (6" casters)	20", 21", 24" (51 cm, 53 cm and 61 cm)
PRODUCT WEIGHT	5" (125 mm) casters: 68 lbs (30.9 kg) 6" (150 mm) casters: 69 lbs (31.3 kg)
FRAME AND HARDWARE MATERIAL	Stainless Steel (Lifetime warranty against corrosion)

All measurements taken with standard molded seat and standard arm pad.

Two-year limited warranty covers the frame against manufacturing defects. A one-year limited warranty covers all other components, subject to normal wear and tear exclusions.

Specifications, and options are subject to change without prior notice.