


XXL DIVISYS 2.0[®]



XXL Divisys bed
with light beech ABÉLIA boards
RAL 1001 PVC hand rail and sidewalls.
Frame color RAL 7035.



BED FOR PATIENTS UP TO 270 KG



Bed on castors with Light Beech skirting with AUZENCE II boards and a Light Beech interior Ral 7030 for surround handrail



Medidom II bed dressing, Light Dragon Tree finish, with long sides, wood skirt dressing and fixed feet.

*Incompatible with rails.



New CA40 box
easy access to the connections.
no vibratory and sound disturbances.

3 basic finishes:
On Castors / On castors with skirting / On fixed feet with skirting.

THE XXL DIVISYS BED ACCORDING TO THE ERGOFIX PRINCIPLE. IT IS BASED ON A PATENTED ASSEMBLY AND DISMANTLING SYSTEM.

PREVENTION OF MUSCULO-SQUELLETAL DISORDERS AND ERGONOMICS:

- ▶ Ergonomic assembly and dismantling.
- ▶ Adjustable height from range from 35 to 81 cm.

COMFORT AND HYGIENE

- ▶ Divisible bed base with welded metal slats.
- ▶ New C40 control box reducing sound and vibratory disturbances.

SAFETY

- ▶ Through rondo screws for assembling the half bed bases.

DESCRIPTION	DIMENSIONS	ADJUSTABLE HEIGHT	BR	LEG REST	LIFTING CAPACITY	CASTORS	STANDARDS	MOTORIZATION
XXL DIVISYS 2 electric functions	120 x 200 cm	35 to 81 cm	70°	12° manual cranks	SWL ⁽¹⁾ 315 kg mattress 30 kg patient 270 kg accessories 15 kg	Ø 100 mm with brake without thread guard	IEC 60601-2-52 compatible with environments 1, 2, 3, 4 and 5	LINAK IP66
	140 x 200 cm							
XXL DIVISYS 3 electric functions	120 x 200 cm			26° thigh sections (thighs) and 16° feet adjustment by manual cranks				
	140 x 200 cm							

⁽¹⁾ Safe working load

- 2 accessory holders (lengthwise)
- 2 IV stand holders (at the head of the bed)
- 2 lifting pole (at the head of the bed)

Handset, electric and mechanical parts: 5-year warranty*

Boards and accessories: 1-year warranty*

*Under normal conditions of use and excluding wearing parts - please see instructions of use

Photos not contractually binding

WINCARE Group - 4, Le Pas du Château - 85670 Saint-Paul-Mont-Penit - France

Tél : +33 (0)2 51 98 55 64 - Fax : +33 (0)2 51 98 59 07 - E-mail : info@medicatatlantic.fr - www.winnicare.fr

In the interest of constantly improving its products, the WINCARE Group reserves the rights to modify the technical specifications of any product without prior notice

Class I medical device in accordance with European Directive 93/42/EEC

