FEATURES



SPECIFICATIONS	Raz-AT600 (Attendant Tilt 600 lbs capacity)
WEIGHT CAPACITY	600 lbs (275 kg)
FRAME WIDTH (seat and back frame)	22" (56 cm)
SEAT WIDTH (commode seat)	22" - 30" (56 cm - 76 cm)
SEAT DEPTH (front of back post to front of seat)	18" (46 cm)
OVERALL WIDTH (with standard armpad)	Seat frame width + 4.5" (11 cm)
OVERALL LENGTH With Z709 Footrests	38" – 39" (96 – 99 cm)
CLEARANCE HEIGHT UNDER SEAT	with pan brackets: 18" - 19" (43 - 48.5 cm) without pan brackets: 18.75" - 19.75" (47.5 - 50 cm)
ARMREST HEIGHT ABOVE SEAT	9.5" (24 cm)
SEAT-TO-FLOOR HEIGHT	21.75" and 22.75" (55 cm and 58 cm)
PRODUCT WEIGHT	5" casters: 82 lbs (37 kg) 6" casters: 84 lbs (38 kg)
FRAME AND HARDWARE MATERIAL	Stainless Steel (Lifetime warranty against corrosion)

All measurements taken with standard molded seat and standard arm pad.

Two-year limited warranty covers the frame against manufacturing defects. A one-year limited warranty covers all other components, subject to normal wear and tear exclusions.

Specifications, and options are subject to change without prior notice.